

KIMA µl TEST

STERILE TUBES FOR CAPILLARY BLOOD COLLECTION



Codes:

- 813005 – Micro test tube **Sterile** K3 EDTA 250 µl
 - 813010 – Micro test tube **Sterile** K3 EDTA 500 µl
 - 8149440 – Micro test tube **Sterile** No Additive 250 – 500 µl
 - 811020 – Micro test tube **Sterile** Clot Activator 500 µl
 - 810020 – Micro test tube **Sterile** Gel+clot activator 800 µl
 - 813805 – Micro test tube **Sterile** KF+Na₂ EDTA 250 µl
 - 812550 – Micro test tube **Sterile** Gel+Lithium Heparin 800 µl
- Box of 50 pieces**

New micro tubes for collection of capillary blood by skin puncture.

The tubes have an easy screw cap closure.

The shape tube's rim have been designed for easy sample collection.

“Micro test” is the unique capillary blood collection tube sterilized by irradiation for maximum safety of patient.



Vacutest Kima S.r.l. | Via dell'industria 2 35020 Arzergrande (PD) Italy
Tel. +39 049 9720624 Fax. +39 049 9720182 www.vacutestkima.it

VACUTEST®
kima