

NEW

URINE CONTAINER REF 90476
CONTENITORE PER URINA COD. 90476

URINE CONTAINER FOR 24 HOURS COLLECTION WITH BUILD-IN STRAW. AMBER COLOR. CAPACITY 3000 ML (GRADUATED MIN. 100 ML – MAX 3000 ML).

CONTENITORE CON SONDA INCORPORATA PER LA RACCOLTA DELLE URINE NELLE 24 ORE. COLORE AMBRATO. CAPACITA' 3000 ML (GRADUAZIONE MIN. 100 ML – MAX 3000 ML).



THE CYLINDRICAL CONTAINER FEATURES A URINE COLLECTION DEVICE FOR A SAFE URINE TRANSFER INTO VACUUM TUBES AND IS PROVIDED WITH AN HANDLE FOR A COMFORTABLE TRANSPORT.

IL CONTENITORE E' DI FORMA CILINDRICA, DOTATO DI DISPOSITIVO PER IL TRASFERIMENTO IN SICUREZZA DELLE URINE NELLE PROVETTE CON VUOTO PREDETERMINATO. IL CONTENITORE E' DOTATO DI MANICO PER UN AGEVOLE TRASPORTO.

THE BUILD-IN STRAW HAVE BEEN DESIGNED TO ALLOW A PROPER TUBE FILLING EVEN WITH 100 ML OF URINE BY TILTING THE CONTAINER AS INDICATED IN THE DRAWING.

LA SONDA INCORPORATA E' STATA PROGETTATA PER PERMETTERE IL RIEMPIMENTO DI UNA PROVETTA ANCHE CON 100 ML DI URINA ATTRAVERSO L'INCLINAZIONE DEL CONTENITORE COME INDICATO DAL DISEGNO.

It. Pat.Pend. PD2009A000373
Registered Design: DM/072971

